


## AN 642

### Oxygen analysis by Flash 2000 Elemental Analyzer

OEA Team

- π Automated, unattended analysis
- π No matrix effect
- π Reliable results
- π Fast analysis



### Introduction

The need of an accurate, fast and automatic method to measure the total Oxygen content in different materials has become of increasing importance in the last years.

The analytical method able to satisfy this need must be such that the experimental data has not to take into account the type of Oxygen bonding in the samples. So, it is necessary to transform in a single chemical species all the Oxygen present in the sample. This species must be able to be detected at high sensitivity due to some of its own physical or chemical characteristics.

In this light the Flash 2000 Analyzer is the most reliable instrument for Oxygen analysis and copes with all requirements of modern laboratories such as accuracy, reproducibility and low cost per analysis.

### Description of the analytical method

The determination of the total Oxygen is based on the Schutz-Unterzaucher technique.

The system operates in pyrolysis mode. Samples placed in silver containers are dropped into the pyrolysis chamber which is maintained at 1060°C and contains nickel coated carbon. The oxygen in the sample combined with carbon forms CO which is then chromatographically separated from other combustion products and sensed by a thermal conductivity detector. Total run time less than 6 minutes.

A complete report is automatically generated by the Eager Xperience data handling software package and displayed at the end of the analytical routine.

### Analytical Conditions

T tube: 1060°C

T oven: 65°C

Helium: Carrier Flow: 130 ml/min  
 Reference Flow: 100 ml/min  
 Sample Delay Time: 12 sec  
 Run time: 300 sec

Standard: CEDFNI, benzoic acid, acetanilide, hexane-ethanol solution (4.5 O%, for gasoline analysis)

Sample weight:  
 Depends of sample nature

Sample pretreatment:  
 Depends of sample nature


## Results

Schema 1 shows the layout of Oxygen configuration.

Table 1 shows a typical analytical test with pure organic compounds. CEDFNI (cyclohexanone-2, 4-dinitrophenylhydrazone, 23 %O) was used as standard, K factor as calibration method and Acetanilide (11.84 %O) was analysed as unknown to check the calibration. The weight of sample was 1-2 mg. The data demonstrate an excellent reproducibility and accuracy.

Table 2 shows Oxygen data of different type of samples. Coal, coke, lignite and compost samples were homogenize using a ball mill before the analysis; while the other samples were analysed without any pre-treatment. Gasoline samples were analyzed by direct liquid injection using a chromatographic syringe. All data were obtained with a good reproducibility and no matrix effect was observed when changing type of sample.

**Analytical layout of Oxygen determination**



**Table 1 – Typical analytical test**

Sample Name	Run	Sample type	Theoretical O%	Experimental O%	Average O%	RSD %
CEDFNI	1	Standard	23.00			
	2		23.00			
	3		23.00			
Acetanilide	4	Unknown	11.84	11.85	11.86	0.4695
	5	Unknown	11.84	11.92		
	6	Unknown	11.84	11.81		

**Table 2 – Oxygen determination of different sample nature**

Sample	Weight (mg)	O %	RSD%
Coal 1	2-3	10.822	0.7951
		10.701	
Coal 2	2-3	3.730	0.5751
		3.723	
		3.690	
Coke 1	1-2	1.044	3.2319
		1.087	
		1.020	
Coke 2	1-2	0.502	1.2012
		0.514	
		0.510	
Lignite 1	2-3	25.202	0.9418
		25.540	
Lignite 2	2-3	21.418	0.7227
		21.200	
Diesel 1	2-3	0.3002	5.0916
		0.2711	
		0.2862	
Diesel 2	2-3	0.2259	6.2487
		0.2560	
		0.2431	
Gasoline 1	2-3 ul	4.0144	0.5865
		4.0477	
Gasoline 2	2-3 ul	0.8497	1.2868
		0.8653	
Gasoline 3	2-3 ul	2.6400	0.1508
		2.6456	

Sample	Weight (mg)	O %	RSD%
Rubber 1	2-3	2.023	0.823
		1.999	
Tyre 1	2-3	4.345	0.254
		4.325	
		4.327	
Wood	3-4	40.223	0.7993
		40.659	
		40.857	
Paper	2-3	42.585	0.4870
		42.918	
		42.969	
Ash	3-4	0.327	1.0730
		0.332	
Carbon fiber 1	2-3	0.2654	1.0382
		0.2679	
		0.2624	
Carbon fiber 2	2-3	7.3561	0.3502
		7.3995	
		7.3536	
Compost	3-4	26.360	1.0165
		26.022	
		26.550	
Stainless steel	20-30	0.0126	1.6129
		0.0124	
		0.0122	
Polymer 1	3-3.5	0.411	1.7032
		0.404	
		0.418	
Polymer 2	3-3.5	0.1595	0.2216
		0.1598	